

Determinação Potenciométrica de Ácido Fosfórico em refrigerante  
06 de novembro de 2013

50 mL de refrigerante + 50 mL de água titulado com NaOH 0,102 mol.L<sup>-1</sup>  
Dados da titulação potenciométrica

NaOH V(mL)	pH
0,0	2,690
1,0	2,910
2,0	3,350
2,5	4,010
3,0	5,750
3,5	6,120
4,0	6,420
4,5	6,710
5,0	6,900
5,5	7,150
6,0	7,450
6,5	7,960
7,0	8,740
7,5	9,220
8,0	9,470
8,5	9,680
9,0	9,800
10,0	10,000
11,0	10,160
12,0	10,250
13,0	10,340
14,0	10,430
15,0	10,490